

Winco Extra Wide Blood Draw Chair | 2580

FACTORY NEW

Winco 2580 - Extra Wide Blood Draw Chair. Winco's Harmony BDC is an extra-wide blood-drawing chair that offers comfort and convenience. The wide seat accommodates bariatric patients, as well as a parent and child. Its design ensures maximum support and ease during blood draws for patients of various sizes.

The chair's adjustable features provide versatility, while the durable construction ensures long-lasting use in medical settings. The padded armrests enhance comfort, and the easy-to-clean upholstery ensures hygiene standards are maintained. Note: Looking for a Blood Chair with different features or weight capacity?

Click the following link to view our catalog of Blood Draw Chairs.



Stock Photo - Options May Vary

FEATURES

- The extra-wide blood-drawing chair offers both comfort and convenience
- Adjustable height, padded arms (option of short arms, L-arms, or a combination of both)
- The wide seat accommodates bariatric patients as well as a parent and a child
- Optional storage cabinet, footrest, and dual armrest padding available

SPECIFICATIONS

- Seat Width: 36.75"
- Weight Capacity: 800 lb
- Seat height: 20"
- Seat Depth: 18.5"
- Height: 35.75"

**Brochure and images are a general overview of the product. Actual product options may vary.*