

# Winco Medical / Recovery / Dialysis Recliner w/Steel Casters | 6531

## FACTORY NEW

Winco 6531 Recliner w/Steel Casters. Outstanding Value and Quality in Our Standard CareCliner. A lower price does not mean lower quality when you buy a Winco CareCliner.

Our standard CareCliner comes with all the superior features you might find on higher-priced clinical chairs. The Winco 6531 CareCliner is our basic workhorse in the clinical recliner line.

Winco CareCliner Features: Winco's 6531 CareCliner has Non-porous blow-molded fold-down side trays mounted on both sides. Each chair in the CareCliner Series features a standard Trendelenburg release with a handle accessible from either side.



*Stock Photo - Options May Vary*

## FEATURES

- Ideal for use in dialysis, blood collection, oncology, respiratory, critical care, and cardiac
- Three-position recliner
- Infinitely adjustable back
- Optional 5" Steel Casters
- Deep recline and Trendelenburg position

## SPECIFICATIONS

- Weight: 94 lbs.
- Height: 49"
- Maximum width: 49"
- Width with Tray Down: 32.25"
- Seat Height: 21"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)