

# Winco Privess Basic 3 Panel Privacy Screen | 3130



Stock Photo - Options May Vary

## FACTORY NEW

The Winco Privess Basic 3-Panel Privacy Screen is a durable steel-framed solution designed for healthcare, disaster relief, and education settings. It features a high-strength powder-coated steel frame with lightweight hinged panels made from flame-retardant 6-mil vinyl. The screen provides portable, cost-effective privacy protection suitable for patient care environments.

The unit includes both casters for mobility and crutch tips for stationary positioning, allowing easy conversion between mobile and fixed configurations. The flat vinyl panel design offers straightforward maintenance and reliable performance in demanding professional settings.

Standard White Vinyl is included, with optional Premium SureChek colored vinyl available for enhanced durability, antimicrobial protection, and stain resistance.

## FEATURES

- 3 Panel Steel Frame
- Comes complete with casters and crutch tips
- High-strength powder-coated steel frame
- Flame retardant 6-mil vinyl
- Easily converts from casters using furnished crutch tips

## SPECIFICATIONS

- Height: 70 in.
- Width: 79 in.
- Weight: 25 lbs.
- Frame Type: Powder-Coated Steel
- Single Panel Width: 27 in.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*