

Champion Privess™ Elite Privacy Screen



Stock Photo - Options May Vary

FACTORY NEW

The Champion Privess Elite Privacy Screen is a professional divider that creates instant personal space in any medical or office setting. It features a durable, powder-coated steel frame for long-lasting use. You can trust this sturdy construction to withstand the demands of a busy clinic.

Consequently, your staff can focus on providing excellent patient care. The smart design offers incredible versatility for your changing floor plan. You may use the smooth casters to move the screen quickly between rooms.

Alternatively, the included crutch tips provide a permanent, stationary base for fixed areas. This flexibility makes it an ideal choice for emergency rooms and triage centers.

FEATURES

- **Three-Panel Design:** Offers a wide 84-inch coverage area to ensure complete patient dignity
- **Steel Frame Construction:** Features a durable, powder-coated frame designed for heavy clinical use
- **Dual Mobility Options:** Includes both twin-wheel casters for transport and crutch-tips for stationary placement
- **Magnetic Locking System:** Keeps panels securely closed and aligned during storage or transport
- **Designer Vinyl Panels:** Uses opaque, commercial-grade vinyl that resists stains and simplifies cleaning

SPECIFICATIONS

- Overall Width (Fully Extended): 84" (213.4 cm)
- Overall Height (With Casters): 70" (177.8 cm)
- Individual Panel Width: 28" (71.1 cm)
- Individual Panel Length (Vinyl): 55" (139.7 cm)
- Base Width: 31" (78.7 cm)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)