

Winco Recovery Couch with Shelf and Contoured Headrest | 8550SH

FACTORY NEW

Winco Recovery Couch w/Shelf & Contoured Headrest: 8550SH. Winco's recovery couches are an excellent choice for almost any medical facility or school nurse's office.

The 8550 Recovery Couch is engineered to provide years of service in any Healthcare setting and is perfect for geriatric or pediatric use. This model includes a paper holder and a cutter and has a wide variety of vinyl colors.

Winco Recovery Couch with Shelf & Contoured Headrest Features Durable foam construction, Sturdy steel frame legs, adjustable leveling feet, and a built-in headrest. The comfortably contoured headrest and 2.5" medium-density foam provide exceptional comfort.



Stock Photo - Options May Vary

FEATURES

- Weight Capacity: 400 lbs
- Height enables the patient to sit or lie comfortably
- Durable foam construction
- Sturdy steel frame legs
- Lower Storage Shelf

SPECIFICATIONS

- Weight: 98 lbs
- Height: 19"
- Width: 28"
- Length: 72"
- Weight Capacity: 400 lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)