

Winco Swing Away Arm Med/Rec/Dia Recliner | 6940

FACTORY NEW

Winco 6940 -Swing Away Arm CareCliner. Design, Comfort, Durability, and Value That Exceed Your Every Expectation. Simple concepts that other manufacturers have yet to grasp.

These concepts set the Winco Elite CareCliner clinical chair series apart from the rest. As a result, elite CareCliners are among the best clinical chairs in the market today.

Winco's 6940 Recliner Features: The wide elastic band on the back frame provides exceptional comfort—standard Trendelenburg release on the back with a handle accessible from either side of the chair. The non-porous blow molded fold-down side tray comes standard on both the right and left sides.



Stock Photo - Options May Vary

FEATURES

- Dual side swing arms for modified entry and ease of cleaning
- 5" nylon total-lock casters
- Built-in headrest Head Rest cover
- Wide full support
- Uni-directional elastic band suspension system integrated into the seat

SPECIFICATIONS

- Chair Weight: 104 lbs.
- Height: 48.5 in.
- Width: 34 in.
- Maximum width: 54 in.
- Width with Tray Down: 34 in.

**Brochure and images are a general overview of the product. Actual product options may vary.*