

Winco Treatment Table with Adjustable Backrest and Cabinet | 8570C1



Stock Photo - Options May Vary

FACTORY NEW

Winco Treatment Table with Adjustable Backrest and Cabinet: 8570C1. For over 60 years, Winco has manufactured furnishings that provide patients with superior comfort. A spirit of innovation and continuous improvement thrives throughout the company and drives its associates to create and introduce “first to market” products to the worldwide healthcare marketplace. Winco’s 8570C1 Treatment Table Features: The Adjustable Backrest is a gas spring operated and uses a ‘one-touch’ adjustment control.

Standard 30” height. Cabinet module and shelf included. High-strength steel frame and legs support up to 400 pounds. The 8570C1’s structure features a powder-coated neutral gray.

2” medium density foam is upholstered in heavy gauge expanded vinyl.

FEATURES

- Standard table
- Laminated storage cabinet
- Lower shelf
- Dove-Gray Powder-Coated Steel Frame
- Paper holder

SPECIFICATIONS

- Weight: 169 lbs.
- Height: 30 in.
- Width: 28 in.
- Length: 72. in.
- Weight Capacity: 400 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)