

# Winco Treatment Table with Adjustable Backrest and Shelf | 8570SH



Stock Photo - Options May Vary

## FACTORY NEW

Winco Treatment Table w/ Adjustable Backrest and Shelf: 8570SH. For over 60 years, Winco has manufactured furnishings that provide patients superior comfort. A spirit of innovation and continuous improvement thrives throughout the company and drives its associates to create and introduce "first to market" products to the worldwide healthcare marketplace.

Winco's Treatment Table w/Adjustable Backrest and Shelf 8570SH Feature: Dove-Gray Powder Coated Steel Frame which supports up to 400 pounds, paper holder, comfortable 2" medium density foam, adjustable leveling feet, and heavy gauge expanded vinyl. The back is gas spring operated and uses a 'one-touch' adjustment control. Standard 30" height.

2" medium density foam is upholstered in heavy gauge expanded vinyl.

## FEATURES

- Adjustable Backrest
- Shelf included
- Gas spring operated
- Uses a 'one-touch' adjustment control
- Standard 30" height

## SPECIFICATIONS

- Weight: 116 lbs.
- Height: 30 in.
- Width: 28 in.
- Length: 72. in.
- Weight Capacity: 400 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)