

Winco Treatment Table with Shelf | 8500SH



Stock Photo - Options May Vary

FACTORY NEW

Winco Treatment Table with Shelf: 8500SH. Customize your Treatment Table! Build onto our basic models using shelves, cabinets, and drawers to meet your specific needs.

The 8500 series are Winco's most versatile and popular treatment tables available. The 8500SH comes with an added shelf for open storage. Choose from a wide variety of beautiful upholstery colors.

Winco 8500SH Features: Dove-Gray Powder Coated Steel Frame, paper holder, cutter, adjustable leveling feet, comfortable 2" medium density foam, and heavy gauge expanded vinyl.

FEATURES

- Added shelf for open storage
- Dove-Gray Powder Coated Steel Frame
- Paper holder and cutter
- Adjustable leveling feet
- Comfortable 2" medium-density foam

SPECIFICATIONS

- Weight: 99 lbs.
- Height: 30 in.
- Width: 28 in.
- Length: 72 in.
- Weight Capacity: 400 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)