

Winco Verō Series 6 Care Cliner™ w/Trend, Swing Arm & 5" Casters | 6Y40

FACTORY NEW

Winco Verō Series 6 Care Cliner™ - 6Y40. Introducing the NEW Verō Care Cliner Collection Designed by Winco, built by you.

That's just smart. Verō Care Cliner 6 Series - Trendelenburg: Winco introduces "Secure-Glide™" Trendelenburg control. Exclusive to the Verō 6 Series, Secure-Glide™ provides gentle patient positioning while preventing clinician strain.

The controlled descent helps eliminate any need, or natural reflex, to strain to control the rapid patient positioning in emergencies. The Verō recliner's modular design helps the clinician choose a configuration that is perfectly suited to their specific needs.



Stock Photo - Options May Vary

FEATURES

- Designed by Winco, built by you
- Provides total flexibility in matching your application with exactly the features you require and the comfort your patients demand
- Promotes patient independence
- Low-effort recline and infinite position backrest
- Supports clinician comfort with bilateral, foot-operated controls for most functions

SPECIFICATIONS

- Weight with Footplate: 114 lbs.
- Overall Height: 48.25 in.
- Overall Width: 30.25 in.
- Seat Depth: 22 in.
- Seat Height (to Floor): 21 in.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)