

Zimmer ATS 1200 Tourniquet System – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

Zimmer ATS 1200 Tourniquet System – Completely Refurbished. This Tourniquet System is a single-port, dual-cuff surgical tourniquet system with microprocessor controls. The Zimmer ATS 1200 is among the safest and most reliable electronic tourniquets.

The ATS 1200 features both audible and visual alarms. In addition, this Zimmer unit features audible and visual alarms that activate when it detects an unsafe condition. It also performs a self-check calibration each time you switch on the system.

Also, this Zimmer 1200 comes completely refurbished with a 6 Months warranty.

FEATURES

- Uses Ambient Air
- Uses Electronic and Microprocessor Regulated Pressure
- Single Port, Dual Cuff System
- Allows for use of IVRA/Bier Block Procedures
- Visual and Audible Alarms

SPECIFICATIONS

- Battery Life: Approx. 45 Minute
- Power Requirement: 120v
- Completely Refurbished

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)