

Zimmer ATS 2000 Tourniquet System – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

Zimmer ATS 2000 Tourniquet System: Refurbished. The Zimmer ATS 2000 is among the safest and most reliable electronic tourniquets. The unit features Digital Displays for Time and Pressure that are large and easy to read. Moreover, fast Inflation and Deflation save time, and Dual Ports make procedures that require bilateral or bier blocks.

Furthermore, Pressure readings are taken at each cuff to ensure accuracy.

The ATS 2000 Also Features: Audible and Visual alarms that will prompt the user of any abnormal conditions with the unit during operation. Furthermore, battery Backup/Operation and a Pole Clamp to mount the unit to a portable stand (portable stand not included).

FEATURES

- Digital Displays for Time and Pressure
- Digital Displays are Large and Easy to Use
- Fast Inflation and Deflation
- Pressure Readings are taken at each cuff
- Audible and Visual Alarms

SPECIFICATIONS

- Dimensions: 12.75"H x 10"W x 10.5"D
- Weight: 19.8 lbs
- Cuff Pressure Range: 50-475 mm Hg, in 5 mm Hg increments
- Time Alarm Set Ranges: 0-240 minutes in 1-minute increments
- Accuracy: +/- 3 mm Hg

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)