

# Zimmer ATS 3000 Tourniquet System – Refurbished



Stock Photo – Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

Zimmer ATS 3000 Tourniquet System – Refurbished. This versatile tourniquet system includes dual ports and cuffs, offering flexibility for a variety of surgical procedures. Equipped with advanced microprocessor controls, it delivers precise pressure management to ensure patient safety and reliable performance.

The fully refurbished ATS 3000 includes a 6-month warranty, offering great value and peace of mind for healthcare providers. In addition, the Zimmer ATS 3000 features audible and visual alarms that activate when an unsafe condition is detected, helping to prevent complications.

It also performs a self-check calibration each time the system powers on, ensuring accurate and consistent operation throughout every procedure.

## FEATURES

- Uses Ambient Air
- Innovative Limb Occlusion Pressure (LOP) Feature
- Uses Electronic and Microprocessor Regulated Pressure
- Dual Port, Dual Cuff System
- Use Dedicated Ports for Supplying and Measuring Pressure

## SPECIFICATIONS

- Battery Life: Approx. 4 Hours
- Power Requirement: 120v
- Completely Refurbished

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)