

Zimmer ATS 4000 Tourniquet System – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

Zimmer ATS 4000 Tourniquet System – Refurbished. This tourniquet system includes dual ports and cuffs, providing versatility for various surgical needs. It features advanced microprocessor controls for precise pressure management, ensuring patient safety and reliable performance.

The ATS 4000 comes fully refurbished and includes a 6-month warranty, offering excellent value and peace of mind. In addition, the Zimmer ATS 4000 Tourniquet System features audible and visual alarms that activate when unsafe conditions are detected. It also performs a self-check calibration each time the system powers on, ensuring accuracy and consistent operation throughout procedures.

This makes the ATS 4000 a trusted choice for surgical teams requiring dependable and efficient tourniquet control.

FEATURES

- Personalized Pressure Technology
- Touchscreen User Interface
- 8.4" Color Display
- Self-check Calibration: automatically checks the accuracy of the machine calibration each time the unit is powered on. This saves setup time.
- Cuff Alert: alerts user of cuff status when attempt is made to set the machine to stand-by.

SPECIFICATIONS

- Height: 13.1 in (33.3 cm)
- Width: 9.3 in (23.6 cm)
- Depth: 10.5 in (26.7 cm)
- weight: 13.2 lbs (6kg)
- Time Alarm Set Range: 1 - 240 minute; 1 minute increments

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)