

Zimmer ATS 750 Tourniquet System - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Zimmer ATS 750 Tourniquet System offers dependable performance and simple operation for surgical teams. This refurbished model provides accurate pressure control and consistent reliability in demanding clinical settings. The lightweight design helps staff move and position the unit with ease.

Additionally, the system supports efficient workflow in operating rooms that require fast setup and clear monitoring. The large LED display delivers clear visibility from multiple angles. The intuitive touchpad lets users adjust time and pressure settings quickly.

As a result, clinicians maintain precise control during procedures. The unit also includes audible and visual alarms, which improve awareness and support safe operation throughout surgery.

FEATURES

- Lightweight
- Accurate
- Easy To Use
- Audible and Visual Alarms
- Battery Backup

SPECIFICATIONS

- Dimensions: 8"H x 7"W x 8"D
- Weight: 6.5 lbs
- Cuff Pressure Range: 50-475 mm Hg, in 5 mm Hg increments
- Time Alarm Set Ranges: 0-240 minutes in 5-minute increments
- Pressure Accuracy: +/- 5 mm Hg

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)