

# Skytron 3600B UltraSlide Surgical Table - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

The Skytron 3600B UltraSlide is a multi-purpose surgical table designed for full-body imaging procedures. It supports a 1,000 lb lift capacity and 800 lb articulated weight limit, making it suitable for a wide range of patient needs. The table features a 23" top slide, wired hand control, and C-arm compatibility for flexible imaging access without patient repositioning.

The table includes hydro-electric brakes and a return-to-level function to enhance workflow and control during procedures. Additional features include accessory side rails compatible with most surgical tools, a removable headrest, and a powered kidney bridge for specialized positioning requirements.

The Skytron 3600B is well-suited for vascular, neuro, urology, and orthopedic procedures. An internal rechargeable battery backup ensures continuous operation during power loss, and the table comes with new pads.

## FEATURES

- Wired Hand Control
- Internal Battery Backup (Rechargeable)
- Return to Level
- Slide Top Movement Capability
- Full Movement Range (Back, Height, +/- Trend, Leg, Tilt)
- Electrically Operated Kidney Bridge
- Excellent C-Arm Access
- 1000 lbs Weight Capacity (Lift), 800 Lbs (articulated)

## SPECIFICATIONS

- Weight Capacity: 1000 lbs Lift, 800 Lbs Articulated
- 23' Slide Top Range
- Table Length (Overall): 82" (with Headrest)
- Table Length: 69" (without Headrest)
- Table Width (Overall): 22"
- Patient Surface Width: 19.5"
- Table Weight: 800 Lbs
- Power Requirements: 120V, 300 Watts

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)