

Steris CMAX Surgical Table – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Steris CMAX Surgical Table – Refurbished. The refurbished Steris CMAX surgical table offers strength, reliability, and modern functionality.

Moreover, with a 1,100 lb weight capacity (in normal position), it supports a wide range of procedures. Additionally, key features include wired hand control, battery operation, hydro-electric brakes, and a sliding tabletop for C-arm compatibility.

It also features Return to Level, North American side rails, and fits most surgical accessories. Because of its modern design and cost-effectiveness, this refurbished surgical table suits many clinical environments.

FEATURES

- C-Arm Compatible Table Top
- Wired Hand Control
- Battery Operation
- Hydraulic Floor Locks
- Return to Level

SPECIFICATIONS

- Weight Capacity: 1100 lbs in Normal Position, 500 lbs in Reverse Orientation.
- Table Length (Overall): 81"
- Table Width (Overall): 21 5/8"
- Tabletop Slide Range: 9" (227 mm) to head, 9" to foot
- Patient Surface Width: 20"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)