

Steris Amsco 3080 RL/SP Surgery Table - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Steris Amsco 3080 RL/SP is a multi-purpose surgical table designed for general operating room use. It features a 400 lb weight capacity and includes wired hand control with battery backup operation. The table is equipped with hydro-electric brakes, C-arm compatibility, and return-to-level functionality.

The table offers comprehensive positioning capabilities including height adjustment, trend, backrest, legrest, and tilt movements, plus manual kidney lift and adjustable headrest. North American side rails accommodate most standard surgical accessories.

This unit has been completely refurbished and includes a new pad set. All mechanical and cosmetic components have been restored to proper working condition.

FEATURES

- C-Arm Compatible Table Top
- Wired Hand Control
- 400 lbs Weight Capacity
- Battery Operation
- Hydraulic Floor Locks
- Return to Level
- North American Side Rails Fit Most Accessories
- Completely Refurbished

SPECIFICATIONS

- Weight Capacity: 400 lbs
- Table Length (Overall): 76"
- Table Width (Overall): 21 5/8"
- Patient Surface Width: 20"
- Pad Thickness: 2"
- Power Requirements: 120V, 4.5 Amps

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)