

Skytron 3500 Slide-Top Surgical Table - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Skytron 3500 Slide-Top Surgical Table is a refurbished operating table with a 500 lb weight capacity, designed for multi-purpose surgical applications including vascular, neuro, urology, and orthopedic procedures. It features wired hand control, a 21-inch top slide capability, hydro-electric brakes, and a return-to-level function for enhanced positioning flexibility.

The table includes North American side rails compatible with most accessories, a removable headrest, and a powered kidney bridge. Its design provides excellent C-Arm access and full-body imaging capability, making it suitable for a wide range of surgical specialties.

This unit has been completely refurbished and comes with new pads. It is equipped with electro-hydraulic floor locks for stability and includes a 6-month warranty.

FEATURES

- Wired Hand Control
- Return to Level
- Slide Top Full Movement Capability
- Electro-Hydraulic Floor Locks
- Excellent C-Arm Access
- 500 lbs Weight Capacity
- 6 Month Warranty

SPECIFICATIONS

- Weight Capacity: 500 lbs
- 21" Slide Top Range
- Table Length (Overall): 79"
- Table Width (Overall): 21.75"
- Patient Surface Width: 19.75"
- Power Requirements: 120V, 300 Watts

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)