

Skytron 6002 Operating Table – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Skytron 6002 is a modern surgical table designed for a variety of procedures including plastic surgery, orthopedics, and general surgery. The table features a 600-lb weight capacity and includes wired hand control for adjusting elevation, leg section, lateral tilt, flex/reflex, Trendelenburg and reverse Trendelenburg positioning, backrest (Fowler), and brake lock/unlock functions. The tabletop rotates 180° and is C-Arm compatible.

The table's headrest is detachable and adjustable via a lever underneath. Four feet automatically engage when using the hand control to secure the table to the floor. The table can be moved on wheels once unlocked via the control button. North American side rails accommodate most standard accessories.

This completely refurbished unit comes with new pads and is ready for immediate use in any surgical setting.

FEATURES

- C-Arm Compatible Table Top
- Wired Hand Control
- 600 lbs Weight Capacity
- 180° Table Top Rotation
- Electro-Hydraulic Floor Locks
- North American Side Rails Fit Most Accessories
- Completely Refurbished
- Comprehensive Warranty

SPECIFICATIONS

- Weight Capacity: 500 lbs
- Table Length (Overall): 76"
- Table Width (Overall): 21¾"
- Patient Surface Width: 20"
- Pad Thickness: 2"
- Power Requirements: 120V

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)