

Skytron 6500 Elite Surgical Table - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

Skytron 6500 Elite Operating Table | Refurbished. The Skytron 6500 Elite is a refurbished surgical table with a 500 lb weight capacity in the standard operating position.

Additionally, features include wired hand control, 180° tabletop rotation, hydro-electric brakes, and C-Arm compatibility. The 6500 also comes with North American Side Rails, which fit most accessories.

This surgical table is ideal for a variety of surgical procedures, including plastic surgery, orthopedics, general surgery, and more. The table utilizes a wired hand control, by which you can adjust the elevation, leg section, lateral tilt, Trendelenburg and reverse Trendelenburg, and backrest.

FEATURES

- C-Arm Compatible Table Top
- Wired Hand Control
- 500 lbs Weight Capacity
- 180° Table Top Rotation
- Electro-Hydraulic Floor Locks

SPECIFICATIONS

- Weight Capacity: 500 lbs
- Table Length (Overall): 76"
- Table Width (Overall): 21¾"
- Patient Surface Width: 20"
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)