

Amsco 2080RC Surgery Table – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Amsco 2080RC Surgery Table | Refurbished. The Amsco 2080RC remote control surgical table is a multi-purpose electric operating table with the convenience of remote control operation. Model 2080 has a 300 lb weight capacity and also features an adjustable headrest.

The Amsco 2080RC Surgery Table has been rigorously and completely refurbished. Amsco 2080RC Features: Back up/down, Leg up/down, Trend/Rev Trend, Table up/down, and Side Tilt. Also included are wired Hand Control (new style), North American Side Rails, which fit most accessories, and Mechanical Floor Locks.

Moreover, includes New Pads and New Labels. IA Option for C-Arm Capabilities Available.

FEATURES

- Wired Hand Control (New Style)
- Mechanical Floor Locks
- North American Side Rails Fit Most Accessories
- 300 lbs Weight Capacity
- IA Option for C-Arm Capabilities Available

SPECIFICATIONS

- Weight Capacity: 300 lbs
- Table Length (Overall): 73"
- Table Width (Overall): 21 5/8"
- Patient Surface Width: 20"
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)