

Amsco 2080M Surgery Table - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Amsco 2080M Surgery Table | Refurbished. The Amsco 2080M manual surgical table is a multi-purpose operating table with the reliability and simplicity of dependable manual operation. 2080 has a 300 lb weight capacity has mechanical floor locks and has gone through our rigorous, complete refurbishing process.

The North American side Rails fit most accessories. Amsco 2080M Also Features: Back up/down, Leg up/down, Trend/Rev Trend, Table up/down, Side Tilt, and Adjustable Head Rest. Includes New Pads and New Labels.

IA Option for C-Arm Capabilities is Available.

FEATURES

- Dependable Manual Operation
- Mechanical Floor Locks
- North American Side Rails Fit Most Accessories
- 300 lbs Weight Capacity
- IA Option for C-Arm Capabilities are available (See Options)

SPECIFICATIONS

- Weight Capacity: 300 lbs
- Table Length (Overall): 73"
- Table Width (Overall): 21 5/8"
- Patient Surface Width: 20"
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)