

# Skytron 6500 HD Hercules Operating Table - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

The Skytron 6500 HD Hercules is a modern surgical table designed for a wide range of procedures including plastic surgery, orthopedics, and general surgery. It features an 850 lbs weight capacity and includes wired hand control for adjusting elevation, leg section, lateral tilt, Trendelenburg, reverse Trendelenburg, and backrest positions.

The table offers 180° tabletop rotation and is C-Arm compatible, making it versatile for various surgical applications. The detachable headrest can be adjusted via a lever underneath. Four feet automatically engage when using the hand control to secure the table in place during procedures.

For mobility and cleaning, the table rolls smoothly on four wheels when the button on the hand control is activated. The unit comes completely refurbished with new pads and North American side rails that accommodate most standard accessories.

## FEATURES

- C-Arm Compatible Table Top
- Wired Hand Control
- 850 lbs Weight Capacity
- 180° Table Top Rotation
- Electro-Hydraulic Floor Locks

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)