

Skytron 3602 UltraSlide Surgical Table - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Skytron 3602 UltraSlide is a multi-purpose operating table featuring a 23" top slide for flexible imaging and positioning. The table maintains a 1,000 lb. lift capacity and 800 lb. articulation capacity regardless of slide position, enabling a wide range of surgical procedures.

The table includes removable leg and back sections for full-body imaging capability, with 9.5" slide-to-head positioning for upper body procedures and 13.5" slide-to-foot positioning for lower body imaging. Accessory side rails accommodate most surgical accessories, and a removable headrest provides additional versatility.

Standard features include a powered kidney bridge, internal rechargeable battery backup, and full movement range including back, height, trend, leg, and tilt adjustments. The table is suitable for vascular, neuro, urology, and orthopedic procedures.

FEATURES

- Internal Battery Backup (Rechargeable)
- Slide Top Movement Capability
- Electrically Operated Kidney Bridge
- Excellent C-Arm Access
- 1000 lbs Weight Capacity (Lift), 800 Lbs (articulated)

SPECIFICATIONS

- Weight Capacity: 1000 lbs Lift, 800 Lbs Articulated
- 23' Slide Top Range
- Table Length (Overall): 82" (with Headrest)
- Patient Surface Width: 19.5"
- Power Requirements: 120V, 300 Watts

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)