

Skytron Elite 6300 Surgical Table – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

Skytron Elite 6300 Surgical Table | Refurbished. Skytron's Elite 6300 surgical table is designed to deliver optimum value, versatility and reliability. It boasts a powerful 1,000 lb. lift and 600 lb.

articulation capacity. Features include 25" Low Height Adjustment, Full Body Imaging Coverage, 210° tabletop rotation, and a 4-point brake system maximizing stability. The 6300 comes with North American Side Rails which fit most accessories.

This surgical table is ideal for a variety of surgical procedures, including: Plastic surgery Orthopedics General surgery and many more.

FEATURES

- Full Body Imaging
- Wired Hand Control
- 1,000 lb. Weight Capacity
- 600 lb. Weight Capacity when articulated.
- 210° Table Top Rotation

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)