

# Skytron 6001 Operating Table – Refurbished



*Stock Photo – Options May Vary*

## FULLY REFURBISHED · CLINICALLY CERTIFIED

The Skytron 6001 is a modern surgical table designed for multi-purpose use across various surgical settings. It features a 500-lb weight capacity in normal position and offers 180° tabletop rotation for optimal patient positioning and surgical access. The table is C-Arm compatible and equipped with North American side rails that accommodate most standard accessories.

This refurbished unit comes complete with a wired hand control, electro-hydraulic floor locks, and new pads. The table has been professionally restored to ensure reliable performance in the operating room.

## FEATURES

- C-Arm Compatible Table Top
- Wired Hand Control
- 500 lbs Weight Capacity
- 180° Table Top Rotation
- Electro-Hydraulic Floor Locks
- North American Side Rails Fit Most Accessories
- Completely Refurbished
- Comprehensive Warranty

## SPECIFICATIONS

- Weight Capacity: 500 lbs
- Table Length (Overall): 75"
- Table Width (Overall): 21 ¾"
- Patient Surface Width: 20"
- Pad Thickness: 2"
- Power Requirements: 120V

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)