

Stryker Trio Surgical Table/Stretcher – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Stryker Trio Surgical Table/Stretcher 1033 functions as both a surgical table and stretcher. It supports the patient from pre-op to recovery, reducing the need for transfers and minimizing injury risks for patients and staff. It features the Quadra Channel™ stabilizing lift system for smooth height adjustments.

The system improves access and maintains patient stability during surgical procedures. Additionally, the Stryker Trio Surgical Table/Stretcher's radiolucent siderails allow fluoroscopy access without moving the patient. This enhances imaging accuracy and streamlines surgical workflows.

Furthermore, four-wheel steel-ring brakes with dual-side activation ensure secure positioning and easy control. Staff can engage brakes from either side.

FEATURES

- Equipped for Use During Pre-op, Procedure, and Post-op
- Quadra Channel™ Stabilizing Lift System
- Siderails Enable Fluoro Clearance
- Four IV Receptacles
- Four-wheel, Steel-ring Brakes with Dual Side Activator

SPECIFICATIONS

- Weight Capacity: 500 lbs
- Table Length (Overall): 80"
- Table Width (Overall): Siderails Up: 33" / Siderails Down: 30"
- Patient Surface Width: 27"
- Siderails: 12" x 49"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)