

Amsco M70 Warming Cabinet - Refurbished



The Amsco M70 Blanket warmer is a cost effective solution for blanket and fluid warming needs in your medical facility.

The M70 features analog controls for the top compartment to adjust temperature, the bottom cabinet stays at constant 120° for blankets. This unit has both visual and audible alarms in case of over-temp detection.

The Amsco M70 warming cabinet comes completely refurbished with a comprehensive warranty and 1 shelf.

- Dual Compartment
- Cost Effective
- Top and Bottom Heating
- Adjustable Temp in Upper Compartment (Dial)
- 120° Temp in Bottom Compartment
- Standard 120V Plug
- Stainless Steel Construction
- Over Temp Alarms (Audible and Visual)
- Includes Blanket Rack and 1 Bottom Compartment Shelf
- Completely Refurbished

Specifications:

- Dims: 26.5" D x 30" W x 74.75" H
- Power: 120V 7.5 Amps
- Standard AC Wall Plug
- Approx Weight: 320 lbs
- Temp Range: 95° - 150°
- Temp Range Upper: 95°-150°
- Temp Lower (Constant) – 120°