

## GE Dinamap Procare 400 Vital Signs Monitor – Refurbished



- Adult or Pediatric Use
- Single Button Operation
- Comfortable for Patients
- Measures NIBP, SpO2, PR, Temp
- Backup Battery
- Ultra Portable
- Rugged
- Recessed Carry Handle
- Clear, Bright Display of Reading
- Alarms
- Lightweight Design
- Built-In Pole Clamp
- Completely Refurbished

The Dinamap Procare 400 Vital Signs Patient Monitor enhances productivity with easy to use features and an intuitive menu structure.

The Procare 400 features NIBP, temp, pulse rate, and SpO2 along with a built in pole mount. With a properly sized BP cuff in place, just press a button to take a BP reading. It's that easy.

The GE Dinamap Procare 400 includes all new accessories (NIBP Hose and Adult Cuff, SpO2 Sensor, Oral Temp Sensor, New Battery and Manual).

This vital signs monitor comes completely refurbished with a comprehensive warranty.

### Specifications:

- Dims: 7.1”H x 8.9”W x 2.4”D
- Weight: 5.7 Lbs
- Battery: Rechargeable
- Battery Run Time: Approx 10-11 Hours, Fully Charged
- Battery Voltage: 6V, 3.3Ah
- Power Requirement: 120v