

Gendron 890 Stretcher - Refurbished



- Trendelenburg/Rev Trendelenburg
- Head Up/Down (Pneumatic Assist)
- Height Up/Down
- Fold Under Siderails
- 4 Wheel Locking
- Steering
- O2 Tank Well
- New Pad
- 6 Months Parts Warranty

Our Completely Refurbished Gendron 890 Stretcher has all of the features you need at a great value! The Easy Access Gendron 890 Stretcher features Trendelenburg and Reverse Trendelenburg, Pneumatic Assist Head Up/Down, easy Height Up/Down, with conveniently located foot pedals on either side.

Additional Gendron 890 Features: Central steering/braking system can be operated from either end/side. Fowler assisted Backrest raises 90° - when in the uprights position, conveniently allows access to the storage area which houses an IV Pole (sold separately) along with security cable and an O2 Tank Well. Full length Siderails fold down and under out of the way to allow easy access. The 890 has large 8" all Surface Carpet Casters - four wheel central casters lock and swivel/steering lock.

Specifications:

- Fowler assisted Backrest: 0- 90°
- 28"x76" patient surface
- Caster Size: 8"
- 4 Wheel Locking Swivel Casters
- 6 Months Parts Warranty