

Hausted 4160 Pediatric Crib - Refurbished



- Includes New Pad
- Hydraulic Height Adjustment
- Manual Backrest
- Plastic Base with Storage
- Includes O2 Tank Well
- Full Length Side Rails
- Central Brake and Steer
- Weight Capacity 325 lbs
- Completely Refurbished

The Hausted 4160 Pediatric Stretcher/Crib is ideal for patients from infant to 6 years old. Three position vertical side rails lock in place where needed and the dual pedestal design provides maximum stability for the patient.

The 4160 is height adjustable and also can be placed into either trendelenburg or reverse trendelenburg. The backrest is manually operated.

The Hausted 4160 Pediatric stretcher/crib comes completely refurbished with a comprehensive warranty and a new pad.

Specifications:

- Minimum Height: 32.5"
- Maximum Height: 39"
- Overall Length: 59.25"
- Overall Width: 29.75"
- Patient Surface: 25.25" x 56"
- Side Rails: 17.75" x 56"
- End Rails: 17.75" x 23"
- Caster Diameter: 8"
- Trendelenburg : 10°
- Fowler Angles: 0-45°
- Weight Capacity: 325 lbs