

Hausted Horizon 462 Airglide Stretcher - Refurbished



- Airglide Tuck Under Side Rails
- Brake and Steer Functions
- \pm Height, \pm Trendelenburg
- Gas Spring Operated Fowler
- 500 lbs Weight Capacity
- Dual Side Foot pedal Controls
- Includes New 4" Pad
- Completely Refurbished

The Hausted Horizon Series AirGlide® Stretcher features include Chrome Plated, Full Length, Tuck Under AirGlide® Patient Safety Rails with Dual Pneumatic Assists, Dual Side Foot Pedal hydraulic Controls for Height and Trendelenburg Adjustment, and Pneumatic Assist Backrest. Brake and Steer Functions can be actuated From Any Corner. The ABS Plastic Base Cover has a storage area and O2 Tank Holder.

The 500 lbs weight capacity can handle large patients with ease.

The Hausted Horizon 462 comes completely refurbished with a comprehensive warranty and a new pad.

Specifications:

- Overall Length: 84"
- Overall Width: 31.5"
- Overall Width (Wide Model): 34.75"
- Patient Surface: 28"W x 76"L
- Patient Surface (Wide Model): 32"W x 77"L
- Height Range: 22.25" – 34.75"
- Weight Capacity: 500 lbs
- Fowler Angle: 0° - 90°
- Trend/Rev Trend: \pm 16°
- Side Rails: 15" x 77"
- Wheel Diameter: 8"
- Pad Thickness: 4"