

Midmark 75L (411) Power Exam Table - Refurbished



- Modern Design
- Hand Control
- Electrically Operated
- +/- Height, +/- Tilt (Trend Only)
- Foot Up/Down and Back Up/Down
- Adjustable Headrest
- Hideaway Adjustable Stirrups
- Pull-Out Cysto Pan
- Table Paper Holder
- Table Paper Straps
- Accessory Side Rails
- Completely Refurbished

The Midmark 75L (411) exam table/chair boasts a very low minimum height of 22.5". This makes it very easy for the patient to get on and off the chair safely and confidently. This power chair comes with a hand control for operating movements, accessory side rail, paper holder, and articulating headrest.

The Midmark Ritter 75L (411) comes with a pull out drain pan at the end of the seat section and pull out heel stirrups. The foot /leg rest section of the 75L is lockable in 4 proctology chair positions.

The Midmark Ritter 75L comes completely refurbished with a comprehensive warranty and new upholstery.

Specifications:

- Minimum Height: 22""
- Maximum Height: 40""
- Maximum Length:
 - Foot Extended: 69.25""
 - Foot & Headrest: 82.75""
- Width: 27""
- Weight Capacity: 325 lbs.
- Paper Roll Holder: 18" x 3.0"
- Hand Control: Height, Tilt, Back & Foot
- Electrical: 115 VAC