

Midmark/Pedigo 500 Crib Stretcher - Refurbished



- Includes New Pad
- One Handed Side Rail Operation: Raise and Lower
- Side rail intermediate stop position
- Raising Backrest
- Trendelenburg / Fowler
- Completely Refurbished

Specifications:

- Top Rail Height: 59 1/2"
- Height with Top: 74"
- Overall Length: 60"
- Overall Width: 30 1/4"
- Spring Height: 32 1/4"
- Caster Diameter: 5"
- Side Rail Intermediate Stop Position: 13 3/4"
- Trendelenburg Angles: - 5,0,5,10 degrees
- Fowler Angles: 0,15,30,45 degrees
- Weight Capacity: 70.0 lbs

The Midmark/Pedigo 500 Pediatric Stretcher/Crib, with its unique oval design, rolls easily on 5" oversize locking casters.

Standard features include: manual fowler, one handed operation lowers and raises side rails with side rail intermediate stop position, three position locking fowler, trendelenburg, reverse trendelenburg, and Raising Backrest.

The Mimark/Pedigo 500 Pediatric stretcher comes completely refurbished with a comprehensive warranty and a new pad.