

Stryker 1069 Eye Surgery Stretcher - Refurbished



- Dual Articulating Head Rest
- Collapsible Side Rails
- Central Brake and Steer
- \pm Height, \pm Trendelenburg
- Gas Spring Assist Backrest
- 500 lbs Weight Capacity
- Includes New 3" Two Piece Pad
- Completely Refurbished

The Stryker 1069 Eye, Head, and Neck Surgical Stretcher offers the greatest head-end clearance in the industry. This stretcher features a dual articulating headrest for precision positioning of the patient during surgery and a gas spring assisted backrest.

This stretcher is ideal for either superior or temporal approach. Positioning controls tuck neatly within the headpiece and allow full articulation. The 1069 also features fold down side rails, four IV receptacles, and four-wheel, steel-ring brakes with foot end brake and steer pedals.

The Stryker 1069 comes completely refurbished with a comprehensive warranty and a new pad.

Specifications:

- Overall Length: 90"
- Overall Width: 31½"
- Patient Surface Width: 26"
- Height Range: 22.5" - 34"
- Weight Capacity: 500 lbs
- Fowler Angle: 0° - 90°
- Trend/Rev Trend: \pm 18°
- Head Adjustment: \pm 4.5°
- Head Piece Tilt: \pm 45°
- Side Rails: 13" x 55"
- Wheel Diameter: 8"
- Pad Thickness: 3"