

Stryker 5050 Stretcher Chair - Refurbished



- Gas Spring Assist Backrest
- Fold Down Side Rails
- Central Brake and Steer
- \pm Height, \pm Trendelenburg
- Built-in O2 Tank Holder
- 400 lbs Weight Capacity
- Includes New 3" Pad
- Completely Refurbished

The Stryker 5050 combines the advantages of a stretcher with that of a recovery recliner. This versatile stretcher chair features a 400 lbs weight capacity, central brake and steer function, O2 tank well and fold down side rails with arm pads. Dual pedestal hydraulics increase stability and the footrest adjusts to various heights to accommodate each individual patient.

Functions include: hydraulic height adjustment, trendelenburg and reverse trendelenburg and a gas spring back rest adjustment.

The Stryker 5050 stretcher chair comes completely refurbished with a comprehensive warranty and a new pad.

Specifications:

- Overall Length: 76"
- Overall Width: 30"
- Patient Surface Width: 24"
- Height Range: 22" - 33.5"
- Weight Capacity: 400 lbs
- Trend/Rev Trend: $\pm 18^\circ$
- Backrest Adjustment: $0^\circ - 90^\circ$
- Side Rails: 10" x 31"
- Wheel Diameter: 6"
- Pad Thickness: 3"