

Steris / Hausted Surgi-Chair ESCEYE-ST - Refurbished



Steris Hausted Surgi-Chair ESCEYE-ST Refurbished. The Hausted Surgi Chair is a multi-position chair and stretcher intended for a patient's treatment, transport, or recovery. With the uniquely designed articulating head section, this unit can be used for specialty procedures that require precision head positioning.

The 500 lbs weight capacity can handle large patients with ease.

The Steris / Hausted SurgiChair ESC-EYE comes completely refurbished with a comprehensive warranty and a new pad.

Features:

- 5" Tente brand casters
- Four-Wheel Brake and Steer Caster System
- Thermoformed base cover
- 24 VDC Backup Battery Operation
- 500 lbs Weight Capacity
- Emergency Backrest release lever
- Includes: New Pad Set, Hand Control, New Batteries, and Manual
- Completely Refurbished

Specifications:

- Overall Length: 75"
- Overall Width: 29.4"
- Patient Surface: 22"W x 75"L
- Height Range: 22.25" – 30.5"
- Weight Capacity: 500 lbs