

## THE 549KL

# DIGITAL PEDIATRIC TRAY SCALE

The 549KL pediatric scale was designed to meet the needs of pediatric practices and visiting clinicians. Its lightweight and compact design ensures convenient mobility. With a maximum capacity of 44 pounds and 20 kilograms, this scale offers providers and patients a range of benefits, including safety, comfort, and flexibility. The scale's curved tray keeps patients calm and secure for accurate readings.

The 549KL offers an easy to clean plastic tray ideal for sterile environments. The user-friendly digital interface includes features such as Tare and Auto-off, making it convenient and easy to use. Its motion detection technology ensures accurate weight measurements even when the patient moves. For easy portability and storage, there is an optional carrying case.



## THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 549KL is your choice for safety, comfort, and reliability.

**120+** LOCAL  
FIELD SPECIALISTS

**30-60 SEC**  
AUTO WEIGHT HOLD

**#1** TRUSTED  
BRAND

**100 YEARS**  
OF SERVICE

**Health o meter®**  
Professional

Health o meter® Professional 549KL

weigh**easier**®

# TECHNICAL SPECIFICATIONS

## Capacity

44 lb / 20 kg

## Resolution

0.5 oz/ 10 g

## Display

0.6", LCD display

## Power Source

9 Volt battery included

## Functions

LB / KG conversion  
 LB / KG lock out  
 Tare  
 Auto zero  
 Auto off  
 Auto hold  
 Motion detection  
 Variable auto-off time  
 Enabling unit view with unit of measure locked  
 Auto-weight hold for 30 to 60 seconds  
 Everlock® - permanent weight lock

## Tray Size (w x d x h)

19 3/8" x 12 1/8" x 2"  
 (492 mm x 308 mm x 51 mm)

## Product Footprint

19 3/8" x 13 1/2" x 4 3/8"  
 (492 mm x 343 mm x 111 mm)

## Product Weight

4 lb / 2 kg

## Portability

64769- optional carrying case  
 2210CART- optional rolling cart

## Warranty

2-year limited manufacturer warranty

# ORDERING INFORMATION

MODEL	DESCRIPTION
549KL	Digital Pediatric Tray Scale, LB & KG
SS-549KL	ScaleSurance® 2 year extended warranty for 549 series*
2210CART	Rolling scale/ medical device cart for 549 series of scales
64769	Carrying case for 549 series of scales



**64769**  
 Carrying Case for 549 Series  
 of Scales



**2210CART**  
 Medical Device Cart for 549  
 Series of Scales

To learn more, call **1-800-815-6615**  
 or visit **homscales.com/549KL**

**Health o meter®**  
 Professional

\*Not available in all countries.