

## Blickman 7921SSD Warming Cabinet -Refurbished



- Single Compartment
- Stainless Steel Construction
- 4 Adjustable Shelves
- PID Digital Controls
- Full Numeric Keypads
- Over Temp Alarms
- Completely Refurbished

The Blickman 7921SSD warming cabinet is a full size warmer that can heat blankets and/or liquids. The unit features a single compartment with four adjustable shelves. The temperature range can be selected between 90° and 160° F. In case of overheating, the unit will alarm, both visually and audible, and shut the heating elements off.

The Blickman 7921 is completely refurbished and comes with a comprehensive warranty.

### Specifications:

- Overall Dimensions:  
30"W x 24" D x 72"H
- Interior Depth: 21 1/2"
- Shelves: 4
- Amps: 8
- Approx Weight: 365 lbs
- Temp Range: 90° - 160°