

Hill-Rom 8050 OB-GYN Stretcher - Refurbished



- Gas Spring Assist Backrest
- Collapsible Side Rails
- Central Brake and Steer
- \pm Height, \pm Trendelenburg
- Built-in O2 Tank Holder
- 500 lbs Weight Capacity
- Includes New 4" Pad
- Completely Refurbished

The Hill-Rom 8050 OB-GYN stretcher offers outstanding patient positioning, quickly and easily. The stretcher comes with integrated foot supports and a sliding top. The foot section retracts out of the way during procedures, leaving the foot rests exposed.

Standard features include: lift assist backrest (fowler), trendelenburg, reverse trendelenburg, central brake and steer, height adjustment, and 500 lbs weight capacity.

The Hill-Rom 8050 OB-GYN stretcher comes completely refurbished with a comprehensive warranty and a new pad.

Specifications:

- Overall Length: 80.5"
- Overall Width: 35.5"
- Patient Surface Width: 30"
- Height Range: 24 $\frac{1}{4}$ " - 38"
- Weight Capacity: 500 lbs
- Pad Thickness: 4"
- Side Rail Length: 37"
- Caster Size: 8"