

# Physio-Control LifePak 20 Defibrillator Biphasic | Pacing | AED - Refurbished



- Biphasic Waveform Technology
- Pacing Function
- AED Function
- Synchronized Cardioversion
- 3 Lead EKG
- Onboard Thermal Printer
- Hands Free Multi-Function Cable
- Built-In Battery and Charger
- Small, Lightweight Design
- Carry Handle
- Many Optional Functions Available
- Completely Refurbished

The Lifepak 20 defibrillator/monitor is the ideal crash cart defibrillator, putting early, effective defibrillation into the hands of BLS responders while transitioning easily to full therapy and monitoring capabilities for the ALS trained code team.

The Lifepak 20 will come with Biphasic waveform technology, non-invasive pacing and AED mode.

Included with the unit: Multi-Function Cable, ECG Cable and Lead wires, Set of Adult Therapy Pads, Printer Paper and New Battery.

Optional accessories and functions include: SpO<sub>2</sub>, Carrying Case and External Paddles.

The Lifepak 20 defibrillator comes completely refurbished with a comprehensive warranty.

## Specifications:

- Max Joules: 360j
- Joule Range: 2-360
- Charge Time: Less than 7 seconds
- Pacing Type: Demand or Non-Demand
- Display Size (Screen): 4.53"W x 3.4"H
- Battery: Rechargeable, Sealed Lead Acid
- Recorder Paper: 50mm, Thermal
- Height (unit only): 21.3 cm (8.4")
- Width (unit only): 26.2 cm (10.3")
- Depth (unit only): 26.2 cm (10.3")
- Weight (unit only): ~5.58 kg (~12.3 lbs)
- Power Requirement: 120v