



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Steris/Amsco Century V116 Sterilizer Refurbished



The Steris Amsco Century V116 sterilizer offered by Venture Medical is a Prevac unit, which is designed for fast, efficient sterilization of surgical instruments and related items. This sterilizer boasts a large 16"x16" chamber size with a depth of 26". The unit is user friendly and has enhanced ease of use with its intuitive touch screen interface screen. The vertical sliding door, which is activated by a foot pedal, makes loading and unloading the unit much easier than with older models that utilized a traditional swinging door.

All refurbished units from Venture Medical come standard with a new 30kW steam generator, unless house steam is available. The Steris Amsco V116 can come in freestanding or recessed mounting configurations and comes with one loading rack and two shelves.

The Steris Amsco V116 sterilizer comes fully refurbished with a comprehensive warranty.

Optional accessories and functions include: Manifold Rack, Loading/Transfer Cart, and Pure Rinse.

- Vertical Sliding Door
- 16" x 16" x 26" Chamber Size
- Touch Screen Display
- Loading Rack and Two Shelves
- Prevac Operation
- Multiple Preprogrammed Cycles
- Self Tests and Alarms
- Help Screen for Troubleshooting
- Steam or Electric Operation
- Single Door
- Freestanding or Recessed Mounting
- Thermal Printer

Specifications:

- Exterior Dims: 35.75"D x 26"W x 74.5"H
- Chamber Size: 16"x 16" x 26"D
- Control Panel Electrical Requirements: 120 Volt, 1-Phase, 9.5 Amps
- Steam Generator Electrical Requirements: Electric Heated – 3 Phase - Call for Details
- Water Feed: 1" NPT, 30-50 psi
- Drain: 1-1/2" ODT
- Operating Temperature: 250° - 270°F
- Recorder Paper: 2.25" W
- UL Listed

VENTUREMEDICAL.COM