

Stryker 1001 Advantage Stretcher - Refurbished



- Gas Spring Assist Backrest
- Glideaway Side Rails
- Central Brake and Steer
- \pm Height, \pm Trendelenburg
- Built-in O2 Tank Holder
- 500 lbs Weight Capacity
- Includes New 4" Pad
- Completely Refurbished

The Stryker 1001 Stretcher Stretchers are perfect for emergency departments, transport, PACU, procedures, or anywhere else a versatile stretcher is needed. The 1001 has a flexible design with 26" patient surface can be configured to best suit your department's needs.

Standard features include: crank backrest (fowler), trendelenburg, reverse trendelenburg, and height adjustment. The 500 lbs weight capacity can handle large patients with ease.

The Stryker 1001 stretcher comes completely refurbished with a comprehensive warranty and a new pad.

Specifications:

- Overall Length: 83"
- Overall Width: 31.5"
- Patient Surface Width: 26"
- Height Range: 22" – 35.5"
- Weight Capacity: 500 lbs
- Trend/Rev Trend: \pm 18°
- Side Rails: 14" x 54"
- Wheel Diameter: 8"
- Pad Thickness: 4"