

Stryker 1080 C-Arm Stretcher - Refurbished



- Full-Length C-Arm Access
- Fold-Down Side Rails
- Central Brake and Steer
- \pm Height, \pm Trendelenburg
- Pneumatic Assist Fowler (x 2)
- End Foot Controls
- 500 lbs Weight Capacity
- Includes New 2" Pad
- Completely Refurbished

The Stryker 1080 Fluoroscopy Stretcher is perfect for patients and procedures that require C-Arm imaging.

The Stryker 1080 features dual pneumatic assist backrests, U-shaped base for ease of access, and a radiolucent litter top.

Additional features include, Trendelenburg and Reverse Trendelenburg along with Central Brake/Steering, Fold Down Side Rails, and end foot controls.

The Stryker 1080 Fluoroscopy stretcher comes completely refurbished with a comprehensive warranty and a new pad.

Specifications:

- Overall Length: 87"
- Overall Width: 30.5"
- Weight Capacity: 500 lbs
- Patient Surface Width: 25.5"
- Fluoro Window: 19.5" – 65.5"
- Height Range: 24" – 33.5"
- Fowler Angle: 0° - 90°
- Trend/Rev Trend: \pm 8°
- Side Rails: 13" x 55"
- Wheel Diameter: 6"
- Pad Thickness: 2"