

## Stryker Renaissance/Advantage Stretcher - Refurbished



Stryker 1731 Pictured\*

The Stryker Renaissance/Advantage Series Stretchers are perfect for emergency departments, transport, PACU, procedures, or anywhere else a versatile stretcher is needed.

Standard features include: lift assist backrest (fowler), trendelenburg, reverse trendelenburg, and height adjustment. The 500 lbs weight capacity can handle large patients with ease.

The Stryker Renaissance/Advantage series stretcher comes completely refurbished with a comprehensive warranty and a new pad.

- Gas Spring Assist Backrest
- Glide-Away or Collapsible Side Rails
- Central Brake and Steer
- $\pm$  Height,  $\pm$  Trendelenburg
- Built-in O2 Tank Holder
- 500 lbs Weight Capacity
- Includes New 4" Pad
- Completely Refurbished

### Specifications\*:

- Overall Length: 83"
- Overall Width: 31.5"
- Patient Surface Width: 26"
- Height Range: 22" - 35½"
- Weight Capacity: 500 lbs
- Trend/Rev Trend:  $\pm$  18°
- Side Rails: 14" x 54"
- Caster Size: 8"
- Pad Thickness: 4"

\*Specifications are approximate.

\*Specifications will vary by model and features selected.