

## Valleylab Force 1C Electrosurgical Unit - Refurbished



- Digital Display
- Bipolar and Monopolar Modes
- Multipurpose ESU
- Compact
- Accurate and Consistent Cutting
- 3 Blended Cut Modes
- REM Fault Indicator
- Accurate and Consistent Cut Performance
- Audible Alarms
- Monopolar Foot-Switch Included
- Optional Bi-Polar Foot-switch
- Optional Cart Available
- Completely Refurbished

The Valleylab Force 1C electrosurgical unit (ESU) is a versatile, multipurpose ESU that is suitable for a wide range of procedures. The unit is compact and a great choice for endoscopic procedures.

The Force 1C features 3 preset blend modes and enhanced safety with the REM Return Electrode Monitoring System. The REM system is vital in helping to prevent electrode burns. The unit also features 200 watts of max cut power, 75 watts in coag, and 50 watts in microbipolar. The easy to read digital display shows the settings and audible alarms alert the user to problems.

Included with the unit: Monopolar foot switch, power cord, operators manual.

The Valleylab Force 1C comes completely refurbished with a comprehensive warranty.

### Specifications:

- Monopolar Cut: 200 Watts
- Blend 1: 175 Watts
- Blend 2: 150 Watts
- Blend 3: 125 Watts
- Monopolar Coag: 75 Watts
- MicroBipolar: 50 Watts
- Height (unit only): 15 cm (6")
- Width (unit only): 28 cm (11")
- Length (unit only): 43 cm (17")
- Weight (unit only): 6.8 kg (15 lbs)
- Power Requirements: 120V, 5A Max