

## Stryker Trio Mobile Surgical Platform/Stretcher - Refurbished



- All-In-One Efficiency
- Manual Operation
- Glide-away Side Rails
- Central Brake and Steer
- $\pm$  Height
- +/- Backrest
- $\pm$  Trendelenburg
- +/- Lateral Tilt
- 500 lbs Weight Capacity
- Includes New 3" Pad
- Completely Refurbished

The Stryker Trio (1033) mobile surgery platform virtually eliminates the need for patient transfers by allowing the patient to stay on the surface for all phases – Pre-Op, Surgery, and Post-Op. The Stryker Trio allows the user to move the unit in the common surgical features, which make it ideal for the ASC and in-office procedure environment. The Trio has manual operations, allowing for simplicity and longevity, as well as a secure braking system, including a steer function for when transporting the patient. mobility, simplicity, and versatility.

Additional features include: standard accessory siderails, glide-away patient siderails, IV pole receptacles, 4-wheel central brake system, 5<sup>th</sup> wheel steering, and low profile base. Optional accessories include drain pan, siderail pads, foot extender, headpiece extender, Eye surgery backrest, and more.

The Stryker 1033 Trio surgical table/stretcher comes completely refurbished with a comprehensive warranty and a new pad.

### Specifications:

- Overall Length: 80"
- Overall Width: 33"
- Patient Surface Width: 27" x 78"
- Weight capacity: 500 lbs
- Height range: 26" - 40"
- Backrest: 1° to 80°
- Lateral Tilt: +/- 15°
- Foot Section: 0° - 50°
- Trend/Rev Trend: +20°/-15°
- Side Rails: 12" x 49"
- Caster Diameter: 6"