

■ SELECT COLLECTION

Augustine

Available: Manual XL



The Augustine is designed with a person-centric approach, to be comfortable and help put the patient at ease. This economical treatment recliner has the clinical features needed for high-use medical environments while engineered to mimic home furnishings with concealed casters and a spacious seat that provides luxurious comfort.



Key Benefits

- Enhanced range of motion for optimum, patient-selectable positioning promotes independence and a relaxing experience.
- Quick Release Seat is easily removed for efficient and thorough cleaning.
- Standard rear central caster lock offers impressive stability and an optional Central Front Brake allows secure braking from the front of the chair, with no obtrusive levers or mechanisms to impede staff or visitors.
- For quality engineering and patient-centric innovation that prioritizes functionality, safety, and comfort, Champion products comply with ISO 13485 and are made in the USA.

Common Applications:

Labor & Delivery Rooms, General Patient Room Seating

Standard Features



Manual Recline



Fixed Arms



3" (7.62 cm)
Dual Casters
With Faux Feet




Central Front
Brake System




500 lbs. (226 kg)
Weight Capacity

Space Requirements


Front View



Top View



Side View



Specifications

Model Number	5001
Manual Recline	Yes
Weight Capacity	500 lbs. / 226 kg
Back Height (from floor ♦ from seat)	45.75" ♦ 28.50" 116.21 cm ♦ 72.39 cm
Overall Width (without tables)	34.50" / 87.63 cm
Overall Width (2 optional tables down ♦ up)	38.50" ♦ 56.50" 97.80 cm ♦ 143.51 cm
Seat Width	26" / 66.04 cm
Seat Height ♦ Seat Depth	19.75" ♦ 20.75" 50.17 cm ♦ 52.71 cm
Arm Height (highest point from floor)	26.75" / 67.90 cm
Overall Upright Depth (with handle)	40" / 101.60 cm
Maximum Recline Depth	76" / 193.04 cm
Dual Zinc Plated Steel Casters (diameter)	3" / 7.62 cm

Available Options




Black Arm
Caps



Heat + Massage
Interior Arm Switch



Headrest Cover
& Pillow Set



Foley Bag
Hook Accessory



IV Pole
Options

(Additional available options are also offered. Talk to a Sales Rep for details.)