

## Luvis M200 LED Surgical Light - New



- Excellent Illumination
- Ergonomic & User-Friendly Design
- Detachable & Autoclavable Main Handle
- No UV Emission & Cool Light
- Excellent Depth of Illumination
- 4,300k Color Temperature
- Intuitive and Easy to Use Controls
- Rotating Central Handle for Focus
- Intensity: 120,000 lux
- Color Rendering Index: 95
- Focusable Light Head
- Lightweight Suspension System
- Various Installation Options
  - Ceiling, Mobile, Dual Ceiling Mount

The Luvis M200 LED Surgical Light is the perfect light for any procedure. The M200 features exceptional quality with energy efficient operation.

The 54 LED'S in each light head produce up to 120,000 lux in a pattern that can be adjusted from 7.8" to 11.8". The Luvis Surgical Light also features a control panel that is synchronized with sub-handle for the user to approach the control panel easily during surgery.

The Luvis M200 surgical light can come in a single, dual or mobile configuration and includes a 2 year mfg. warranty.

### Specifications:

- LED's: 54
- CRI: 95
- Color Temperature: 4,300k
- Head Size: 19.6"
- Focus Size: 7.8" (Circle)/11.8" (Circle)
- Focal Distance: 1M
- Intensity: Max. 120,000lux
- Input Power: 100-240Vac, 60Hz
- Power Consumption: 85W
- LED Life Span: 50,000hr